


# Comfort Series 300



Medispecs Comfort Series 300 Incorporates 3 electronic motors that work together to control height, back rest, knee-break and calf section adjustment.



## Comfort Series 300 HB03-S

Medispecs Comfort Series 300 has been designed to maximise patient positioning, comfort and accessibility. The three motor HB03-S is our entry level bed which offers amazing quality and functionality at an affordable price. Ideal for use in clinics, surgical, medical and general wards

### Key Features:

- Electrically operated height, back rest, knee-break and calf section adjustment.
- Trendelenburg and reverse Trendelenburg is controlled manually.
- Dual sided manual CPR.
- Battery Back Up.
- Antimicrobial, waterproof mattress.
- 4 hook Intravenous pole.
- Angle Indicator for back rest and thigh section.
- Side rails are removable for easy cleaning and manufactured from vacuum formed high impact plastic.
- Illuminated patient hand remote control .
- Safe working load of 250kg.
- Variable height of 44cm to 76cm.
- Anti Static directional castors with linked braking system.
- Customisable trim plates. (Colours/Logos/Patterns)

### Technical Specifications

Width	: 104 cm
Length	: 220 cm
Width (Lying Surface)	: 90 cm
Length (Lying Surface)	: 196 cm
Back Segment Length	: 68 cm
Back Segment Angle	: 0-70°
Height Adjustment	: 44-76 cm
Foot Segment Angle	: 0-45°
Trendelenburg	: 12° (Manual)
Safe Maximum Load	: 250 kg
Shipment Dimension	: 220x113x90 cm



### Standard Accessories

- Infusion Pole
- Drainage Pole
- Mattress

### Optional Accessory

- X ray Compartment



### Technical Support



### Quality

