


Comfort Series 400



Medispecs Comfort Series 400 Incorporates 4 electronic motors that work together to control height, tilt, back rest, knee-break and calf section adjustment.



Comfort Series 400 HB04-S

Medispecs Comfort Series 400 has been designed to maximise patient positioning, comfort and accessibility. The four motor HB04-S is ideal for acute care patients as well as medical and surgical wards.

Key Features:

- Electrically operated height, tilt, back rest, knee-break and calf section adjustment with electrically operated Trendelenburg and reverse Trendelenburg.
- Electronic Auto CPR.
- Battery Back Up.
- Antimicrobial, waterproof mattress.
- 4 hook Intravenous pole.
- Angle Indicator for back rest and thigh section.
- Side rails are removable for easy cleaning and manufactured from vacuum formed high impact plastic.
- Illuminated patient hand remote control .
- Calf section can be lowered and raised manually by Carer.
- Safe working load of 250kg.
- Variable height of 37cm to 77cm.
- Anti Static directional castors with independent braking system.
- Customisable trim plates. (Colours/Logos/Patterns)

Technical Specifications

Width	: 104 cm
Length	: 220 cm
Width (Lying Surface)	: 90
Length (Lying Surface)	: 196
Back Segment Length	: 68 cm
Back Segment Angle	: 0-70°
Height Adjustment	: 37-77 cm
Foot Segment Angle	: 0-45°
Trendelenburg	: 17°
Reverse Trendelenburg	: 17
Safe Maximum Load	: 250 kg
Shipment Dimension	: 220x113x90 cm



Standard Accessories

- Infusion Pole
- Drainage Pole
- Mattress
- Battery Back Up

Optional Accessories

- Manual CPR
- Linked Braking system
- X ray Compartment

Technical Support



Quality

