


# Delivery/Gynaecology Bed



Medispecs Delivery Bed Incorporates 3 electronic motors that work together to control height, tilt and back rest adjustment.



## Delivery/Examination Bed HG02-S


Medispecs multi purpose birthing and maternity bed combines inspection, delivery and post natal functions in a single unit. The innovation in design allows this bed to extend from a Gynaecology table to a full sized bed with the aim of a smooth labour process for both the patient and the caregiver.

### Key Features:

- Electrically operated height, tilt, back rest, Trendelenburg and reverse Trendelenburg.
- Detachable foot section allows the bed to be used as a standard bed as well as a delivery bed.
- Battery Back Up.
- Antimicrobial, seamless waterproof mattress.
- Stainless steel waste bin.
- 4 hook Intravenous pole.
- Angle Indicator for back rest.
- Side rails are removable for easy cleaning and manufactured from vacuum formed high impact plastic.
- Illuminated patient hand remote control .
- Safe working load of 200kg.
- Variable height of 66cm to 96cm.
- Anti Static directional castors with independent braking system.

### Technical Specifications

Width	: 105,5 cm
Length	: 201 cm
Back Segment Length	: 85 cm
Back Segment Angle	: 0-70°
Foot Segment Length	: 58 cm
Height Adjustable	: 66,5-96,5 cm
Trendelenburg	: 14°
Reverse Trendelenburg	: 14°
Safe Maximum Load	: 200 kg
Net Weight Gross	: 182 kg
Weight Shipment	: 240 kg
Shipment Dimension	: 220 x 113 x 90cm
Material car	: 66 x 112 x 99cm



### Standard Accessories

- IV Poles
- Knee support
- Heel Support
- Waste Container
- Anesthesia Frame
- Mattress
- Material cart

### Technical Support



### Quality

