

# Dialysis Couch



Medispecs Dialysis Couch Incorporates 4 electronic motors that work together to control height, tilt, back rest, knee-break and foot section adjustment.



## Dialysis Couch DB01-L

Medispecs Dialysis couch has been specifically designed to maximise comfort and accessibility for patients receiving Dialysis treatment by providing for both seated and reclined positions.

### Key Features:

- Electrically operated height, back rest, knee-break and calf section adjustment with electrically operated Trendelenburg.
- Automatic CPR/Emergency function.
- Battery Back Up.
- Antimicrobial, waterproof mattress.
- 4 hook Intravenous pole.
- Side rails are removable for easy cleaning and manufactured from vacuum formed high impact plastic.
- Patient/Carer hand remote control .
- Safe working load of 250kg.
- Variable height of 55cm to 87cm.
- Anti Static directional castors with independent braking system.
- Customisable trim plates. (Colours/Logos/Patterns)

### Technical Specifications

Width	: 93 cm
Length	: 215 cm
Back Segment Length	: 70 cm
Back Segment Angle	: 0-75°
Height Adjustment	: 55-86,5 cm
Foot Segment Angle	: 0-35°
Trendelenburg	: 15°
Safe Maximum Load	: 250 kg
Shipment Dimension	: 220x113x90 cm



### Standard Accessories

- Infusion Pole
- Drainage Pole
- Mattress

### Technical Support



### Quality

