



MONEPLUS GENERAL PROCEDURE CHAIR BA11-L



Medispecs Moneplus General procedure and Transfer chair provides side rails and independent locking anti static castors for safety and ease of transfer between procedures. Fully motorised with 4 electric motors and Trendelenburg, this transfer procedure chair is versatile and effective.

Key Features

- Electrically operated height, tilt, back rest, and knee-break with electrically operated Trendelenburg.
- Antimicrobial viscoelastic foam is used for the padding. This seamless padding moulds itself to the body providing additional comfort and support for long periods. Seamless padding reduces the risk of cross contamination by eliminating sewing and creating a smooth clear surface area which is easy to clean and disinfect.
- Detachable armrests and high impact side rails.
- Hand remote control
- 355 Degree rotation.
- The detachable foot rest and head rest is fully adjustable to accommodate the height of the patient.
- Safe working load of 230kg.
- Variable height of 48cm to 78cm.
- Anti Static directional castors with independent braking system.
- This multi purpose chair can be used in for all general procedures and transfers.



Technical Specifications

Width	: 92 cm
Length	: 183 cm
Back Segment Length	: 68 cm
Back Segment Angle	: 0-80°
Min-Max	: 48-78 cm
Seat Length	: 53,7 cm
Foot Segment Length	: 22,5 cm
Foot Cramp Support	: 13 cm
Foot Segment Angle	: 0-90°
Trendelenburg	: 0-24°
Safe Maximum Load	: 230 kg
Net Weight Gross	: 128 kg
Weight Shipment	: 158 kg
Dimension	: 161x91x88 cm

Standard Accessories

- Plastic Side Guards
- Anti Static Castors
- Head Cushion

Optional Accessories

- Battery Back Up
- Soft Cushion

Color Options



Cushion Accessories



TechnicalSupport



Quality

