



## MONESTAR PODIATRY CHAIR BA04-L



Medispecs Monestar Podiatry Chair comes standard with individually powered legs, comfortable seating and 355 degree rotation. Its design is specific to the industry and suits both patient and caregiver in terms of functionality and comfort. Fully motorised by 4 electric motors with Trendelenburg

### Key Features

- Electrically operated height, tilt, back rest, and individually operated knee-break with electrically operated Trendelenburg.
- Antimicrobial viscoelastic foam is used for the padding. This seamless padding moulds itself to the body providing additional comfort and support for long periods. Seamless padding reduces the risk of cross contamination by eliminating sewing and creating a smooth clear surface area which is easy to clean and disinfect.
- Adjustable armrests.
- Hand remote control
- 355 Degree rotation.
- The detachable foot rests and head rest is fully adjustable to accommodate the height of the patient.
- Safe working load of 230kg.
- Variable height of 48cm to 78cm.
- Anti Static directional castors with independent braking system can be added as an optional extra.
- This Chair can be used for a broad spectrum of foot care procedures.



**Technical Specifications**

Width	: 100 cm
Length	: 185 cm
Back Segment Length	: 68 cm
Back Segment Angle	: 0-80°
Min-Max	: 58,5-90,5 cm
Seat Length	: 51 cm
Foot Segment Length	: 23,5 cm
Foot Cramp Support	: 14 cm
Foot Segment Angle	: 0-90°
Safe Maximum Load	: 230 kg
Net Weight	: 104 kg
Gross Weight	: 134 kg
Shipment Dimension	: 161x91x88 cm

**Standard Accessories**

- Head Support
- Foot Cramp Support

**Optional Accessories**

- Battery Back Up
- Soft Cushion
- Towel-Holder

**Color Options**



**Cushion Accessories**



**Technical Support**



**Quality**

