



MONESTAR PROCEDURE CHAIR BA09-L



Medispecs Monestar Procedure Chair is cost effective and Customizable to specific procedural needs. Fully motorised by 4 Electric motors with Trendelenburg and 355 degree rotation, this chair is versatile and practical for every caregiver and patients needs.

Key Features

- Electrically operated height, tilt, back rest, and knee-break with electrically operated Trendelenburg.
- Antimicrobial viscoelastic foam is used for the padding. This seamless padding moulds itself to the body providing additional comfort and support for long periods. Seamless padding reduces the risk of cross contamination by eliminating sewing and creating a smooth clear surface area which is easy to clean and disinfect.
- Detachable armrests.
- Hand remote control
- 355 Degree rotation.
- The detachable foot rest and head rest is fully adjustable to accommodate the height of the patient.
- Safe working load of 230kg.
- Variable height of 48cm to 78cm.
- Anti Static directional castors with independent braking system can be added as an optional extra.
- This multi purpose chair can be used in numerous disciplines and procedures such as Dialysis, Chemotherapy, eye operations, Oncology and Skin care procedures.




Technical Specifications

| | |
|---------------------|-----------------|
| Width | : 100 cm |
| Length | : 185 cm |
| Back Segment Length | : 68 cm |
| Back Segment Angle | : 0-80° |
| Min-Max | : 58.5-90.5 cm |
| Seat Length | : 51 cm |
| Foot Segment Length | : 23,5 cm |
| Foot Cramp Support | : 14 cm |
| Foot Segment Angle | : 0-90° |
| Trendelenburg | : 0-24° |
| Safe Maximum Load | : 230 kg |
| Net Weight Gross | : 128 kg |
| Weight Shipment | : 158 kg |
| Dimension | : 155x106x90 cm |

Standard Accessories

Color Options



Cushion Accessories



Technical Support



Quality

