



MONESTAR ECO RECLINER CHAIR DA05-E



Medispecs Monestar Eco Recliner chair is affordable, electrically powered and compact. Offering two motors, a hand control, multi positioning arm rests and four locking anti static castors, this chair offers exceptional functionality and comfort. Its compact, low position allows for ease of access for wheel chair patients and transfers

Key Features

- Electrically operated back rest, knee-break with manual adjustable foot section.
- Polyurethane foam is used for the padding. This seamless padding moulds itself to the body providing additional comfort and support for long periods. Seamless padding reduces the risk of cross contamination by eliminating sewing and creating a smooth clear surface area which is easy to clean and disinfect.
- 4 hook Intravenous pole
- Adjustable armrest with left/right rotation
- Hand remote control.
- The foot rest and head rest is fully adjustable to accommodate the height of the patient.
- Safe working load of 170kg.
- Fixed height of 57cm.
- Anti Static directional castors with independent braking system.
- This chair is mainly used for blood transfusion/donation but can also be used for dialysis procedures.



Technical Specifications

Width	: 86 cm
Length	: 178,5 cm
Back Segment Length	: 59.5 cm
Back Segment Angle	: 0-75°
Min-Max	: 57 cm
Seat Length	: 49 cm
Foot Segment Length	: 25.5 cm
Foot Cramp Support	: 15 cm
Foot Segment Angle	: 0-75°
Safe Maximum Load	: 170 kg
Net Weight	: 56 kg
Gross Weight	: 86 kg
Shipment Dimension	: 161x91x88 cm

Standard Accessories

- Intermediate Triangle Pillow

Optional Accessories

- Infusion Pole
- Soft Cushion
- Towel-Holder
- Viscoelastic foam

Color Options



Optional Accessory



Technical Support



Quality

