


Comfort Paediatric Bed



Medispecs Paediatric Bed Incorporates 3 electronic motors that work together to control height, tilt and back rest adjustment.



Comfort Paediatric Bed HB08-S

Medispecs Paediatric bed is perfect for children up to 5ft. in length. This Paediatric bed offers 3 electric motors and an extremely low position of 39cm for special needs and at risk patients. With easy to use side rails, the positioning, accessibility and comfort are unmatched in this category.

Key Features:

- Electrically operated height, tilt and back rest adjustment with electrically operated Trendelenburg and reverse Trendelenburg.
- Electronic Auto CPR.
- Battery Back Up.
- Antimicrobial, waterproof mattress.
- 4 hook Intravenous pole.
- Tilt angle Indicator.
- Side rails are removable for easy cleaning and manufactured from vacuum formed high impact plastic.
- Patient/Carer hand remote control .
- Safe working load of 170kg.
- Variable height of 39cm to 75cm.
- Anti Static directional castors with independent braking system.
- Customisable trim plates. (Colours/Logos/Patterns)

Technical Specifications

Width	: 85 cm
Length	: 166 cm
Width (Lying Surface)	: 75 cm
Length (Lying Surface)	: 155,5 cm
Back Segment Length	: 66,5 cm
Back Segment Angle	: 0-75°
Height Adjustment	: 39-75 cm
Trendelenburg	: -16 / +16°
Safe Maximum Load	: 170 kg
Net Weight Gross	: 158 kg
Weight Shipment	: 208 kg
Dimension	: 220x113x90 cm

Standard Accessories

- Infusion Pole
- Mattress

Optional Accessory



Technical Support



Quality

